

# **Music Box Mechanism**

P7-7330

# **INSTRUCTIONAL GUIDE**

## Contents

- Music Box Mechanism
- Instructional Guide



## Background

Sound travels as a longitudinal wave through matter. Sound waves are a result of mechanical energy from a vibrating mass changing the pressure of the medium the sound is traveling through. By placing a vibrating object against a larger resonant object such as a drum head, chalk board, or thin board, the sound is amplified. The larger surface area of the resonant object creates more sound waves than the vibrating object alone.

## Instructions

## Before Use:

- 1. Unscrew the small screws on top of the plastic case to separate the top and bottom of the case
- 2. Pop the metal mechanism out and attach the handle to the winding mechanism.

## **Demonstrations:**

- 1. Hold the mechanism in your hand (or, better yet, hang it from a string) and wind it. Ask students to listen and describe the loudness of the sound.
- 2. Amplify the sound with one of the methods listed below.
- 3. Ask students to predict which amplification method will produce the loudest sound.

## To amplify the sound, hold the music box mechanism against:

- Chalkboard
- Desk
- Window
- Piece of paper
- Piece of paper rolled into a cone
- Piece of cover stock rolled into a cone

- Drum
- Guitar body
- Piano cabinet
- Student's temple
- Student's elbow, as they press an index finger against the bone near their ear

# **Related Products**

Sympathetic Tuning Fork Set (P7-6000) Set includes two C Note 256 Hz tuning forks, one fixed and one adjustable, each mounted on a wooden resonance box and a mallet.

**Loudspeaker Kit (P7-7800)** This DIY speaker kit is all you need to construct 10 loudspeakers so that students can get real, hands-on experience of how the device works.

**Piezo Buzzer in Vacuum Chamber (97-6600)** A modern replacement for the traditional Bell Jar and Ringer. An acrylic chamber houses a battery-operated, high output piezo sounder. When the end plates are attached and the unit connected to a vacuum pump no sound can be heard but the sound returns when the air is allowed back in.